**PULPDENT®**

# Multi-Cal™

## Calcium Hydroxide Preparation

**Caution:** Federal law restricts this device to sale by or on the order of a dentist.

Made in USA

**For all clinical applications where calcium hydroxide is indicated**

Store in Cool, Dry Place

### • Bactericidal • Bacteriostatic • pH > 12

Pulpdent Multi-Cal is a smooth, creamy calcium hydroxide preparation in a radiopaque hydroxyethylcellulose base and is recommended for all clinical applications where calcium hydroxide is indicated. It is a dressing, not a hard base, and it provides all the same clinical results as Pulpdent Cavity Liner, Pulpdent Paste and TempCanal.

**Indications:** Cavity Lining, Indirect Pulp Capping, Direct Pulp Capping, Pulpal Curettage, Pulpotomy, Dentin Bridge Formation, Temporary and Intermediate Canal Treatment, Periapical Cysts, Abscesses and Lesions, Root Resorption and Apexification

**Note:** To extend shelf life of Multi-Cal, place a drop of water in the cap and recap syringe promptly after each use.

#### CAVITY LINING

Multi-Cal, with pH > 12, neutralizes acids and chemical irritants found in dental filling materials and etch gels. It desensitizes the dentin by stimulating the formation of sclerotic and secondary dentin and protects the pulp against thermal and chemical shock and dehydration.

#### Directions for Cavity Lining and Indirect Pulp Capping:

1. Apply a very thin film to the dentin.
2. Dry with a slow flow of air until material appears chalky white (30 seconds).
3. Remove from the enamel margins.
4. When using glass ionomer and dentin bonding materials, place Multi-Cal in the deepest part of the cavity, leaving the facing walls available for dentin bonding.

#### VITAL PULP TREATMENT

Multi-Cal is applied directly on the pulp in all pulp capping, pulpotomy and pulpal curettage procedures. It stimulates the pulp to form a new dentin bridge which is visible clinically after 2 - 3 months.

#### Directions for Direct Pulp Capping:

1. Apply a very thin film of Multi-Cal on all exposed dentin and dry with a slow flow of air.
2. Place a bead of Multi-Cal 3mm thick directly over the exposed pulp and dry with a slow flow of air until a crust appears. Multi-Cal does not set hard.
3. Gently tamp the paste with a pellet of cotton to make sure it is in contact with the pulp.
4. Place a thin mix of hard setting cement over Multi-Cal.
5. After cement hardens, place final restoration.

**Note:** For pulpotomy and pulpal curettage procedures, pack pulp gently with sterile cotton until blood clots, then proceed as above. A local anesthetic may be used.

#### ROOT CANAL TREATMENT

Use routinely as a bactericidal, bacteriostatic canal dressing between off visits and as a dressing for cysts, abscesses, lesions, etc. Also use for apexification and other hard tissue formation.

#### Directions for Root Canal Therapy:

1. Prepare and cleanse canal in conventional manner.
2. Fill canal with Multi-Cal.
3. Seal access opening with your preferred cement.
4. Multi-Cal is easily removed with file and irrigation.
5. In cases requiring extended root canal therapy, change dressing every month.

#### CLINICAL SITUATIONS WHERE CALCIUM HYDROXIDE CAN BE USED

1. exudation control
2. large periapical lesions
3. antibacterial dressing
4. temporary root filling
5. apical inflammatory resorption
6. inflammatory resorption following trauma
7. apical internal resorption
8. internal-external root resorption

9. perforations
  10. transverse root fractures
  11. incompletely developed pulpless teeth
- (Heithersay, G.S. Calcium Hydroxide in the Treatment of Pulpless Teeth with Associated Pathology, J. Brit Endo Society, 8[2]: 74-93 July 1975)

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***PULPDENT Corporation***